

1. Apply ageLOC Body Shaping Gel to the area to be immediately treated.
2. Turn on the ageLOC Galvanic Body Spa
 - a. Press the on/off button once. The ageLOC rings will appear. Treatment time will display the 5 min.
3. Perform a treatment with ageLOC Body Shaping Gel
 - a. First thoroughly moisten your fingers and palm with water or NaPCA Moisture Mist. Hold the instrument, making sure your fingers are touching the chrome activators. Touch the conductive surface to your skin where the ageLOC Body Shaping Gel is applied.
4. If you lose contact with the chrome activators at any time, simply adjust your hold on the device and resume treatment. Make sure that all contact points are thoroughly moistened.

Recommendations

- For complete usage instructions and precautions, please refer to the ageLOC Galvanic Body Spa Owner's Manual. Be sure to review all directions before using the Body Spa and Body Shaping Gel.
-